


DIESEL ENGINE MODEL

N98-OSK15855




This product is a tabletop model of a diesel engine. It has a transparent body, which allows you to observe the movement of the piston, crank, valve, and rocket arm. Combustion is illustrated by a lamp at the end of each cycle.

WEIGHT AND DIMENSIONS

22cm x 18cm x 34cm (L x W x H)

Net Weight: 3 kg

Gross Weight: 5 kg

OPERATES
ELECTRICALLY

LOOK US UP ONLINE

Version2.1.JL091914