

ATOMA MICROSCOPE

N99-D21-4151

\$495



The versatile ATOMA microscope is specifically designed for advanced high school students and college field research. The heavy-duty construction permits indoor and outdoor observation to occur. This versatile portable model is unique in design and contains a built-in rechargeable battery. When fully charged, the ATOMA can be used up to 1,000 hours. When time is of the essence, a simple 10 minute charge (quick charge) permits 1 hour of usage time. The built-in rechargeable battery has a life expectancy of over 500 charges.



Service Receptacle Connection

ACCESSORIES

- Dust Cover
- Removable AC Cord



LED Rapid Charge Indicator

SPECIFICATIONS

Magnification: 40X-300X,
40X-400X

Objective: 4X, 10X, S40X
Eyepiece: WF10X

Tube: Length 139mm; 45 degrees
inclined and 360 degrees Rotary
Revolver (nose piece): 3 pieces

Stage: M-Stage

Functions: Vertically Movable

Lighting: LED power supply with dimmer and charging lamp

Accessories: Removable AC cord, dust cover

FEATURES

- White LED light for clear observation
- Stage turns 360 degrees
- M-stage moves from front to back, and from side to side
- Easy-focusing function for precise interpretation of results
- Easy-holder function for quick and simple adjustments
- Easy-grip function for portability
- Spot color coding for easy adjustment