





Speed, Performance, Value - Perfect Balance

- Stainless steel pan
- Clear backlit LCD display
- Auto power off
- Mains or battery operated
- Mains adapter supplied as standard
- Several weighing units
- Draft shield supplied as standard
- Integral security bracket
- Adjustable feet
- Weighing



Technical Specifications

Maximum capacity	600g
Readability	0.1g
Repeatability (S.D.)	0.1g
Linearity (±)	0.1g
Pan size	5.1" Ø / 130mm Ø

General Specifications

Model	AQT 600
Resolution	1:6000
Units of measurement	Gram (g), Carat (ct), Ounce (Oz), Pennyweight (dWt)
Typical stabilisation time	2 Sec
Calibration	Push button using external mass
Display	6 digit LCD, 15mm digits with backlight
Power supply	12VAC @ 150mA adapter or 6 x AA batteries
Operating temperature	32°F - 104°F
Other features	Low battery, stable and zero indicators
Housing	ABS plastic
Draught shield dimensions	7.1 x 6.2 x 3.2"
Overall dimensions (wxdxh)	7.7x9.4x2.8" / 195x240x70mm

Key Accessories

- In-use Cover
- Hard Carry Case with Lock
- Adam Lock and Cable



	NADA Scientific
	39 Butternut Street
	P.O. Box 1336
	Champlain, NY 12919
	E-mail: support@nadascientific.com
Tel.: 1-800-799-6232, Fax: 518-297-3524	



Adam Equipment - offices worldwide

www.adamequipment.com