


TURBOJET TURBINE ENGINE MODEL, HIGH PRESSURE COMPRESSOR

N98-ND9210




This super model of a modern two-wave turbine shows in detail the construction and the operating system of this type of motor. The section cut housing allows for a clear view of the high-pressure compressor, high-pressure turbine, combustion chambers with the injection nozzles, and starting plugs. The turbines can be set in motion with the built-in electric motor.

WEIGHT AND DIMENSIONS

60cm x 40cm x 40cm (L x W x H)

Net Weight: 10 kg

Gross Weight: 20 kg

OPERATES
ELECTRICALLY

MOUNTED ON BASE

Version2.1.JL030716



LOOK US UP ONLINE